



FREE Shipping **OxygenPortableUSA**

Your one stop shop for Portable Oxygen Equipment, Service and Accessories

Live Customer Support

9am to 6pm PST

1.800.339.1558

OXLIFE INDEPENDENCE


CONCEPTS



THE PIECE THAT COMPLETES YOUR NEEDS.

The Oxlife Independence by O2 Concepts is a portable oxygen concentrator designed with both patients and providers in mind. Boasting both continuous flow up to 3L/min and pulse settings up to 6 (96 mL), the Oxlife Independence gives patients total freedom with 24/7 oxygen on the go while offering providers a one solution approach to maximize profits in today's turbulent reimbursement climate.

-  Patent-pending ESA technology enables the Oxlife Independence to have the longest battery life in its class

- Only POC in class that can be used on all settings while plugged in car (DC) mode



- Integrated DC power supply -no external brick

- Custom split battery configuration allows for "hot swapping" of battery. No loss of oxygen when switching battery

- Integrated cart kit design offers lighter overall weight

- Integrated cart with 6" wheels offers "best in class" mobility



- Convenient access to preventative maintenance door for fast and easy pump filter change

- Medicare HCPCS codes E1390/E1392

- FAA Approved



PULSE FLOW SETTINGS



- Pulse flow settings in .5 increments for better patient titration
- Large intuitive LCD screen eliminates patient confusion

SETTING	DOSAGE (ML)
.5	8
1	16
1.5	24
2	32
2.5	40
3	48
3.5	56
4	64
4.5	72
5	80
5.5	88
6	96

PRODUCT SPECIFICATIONS*

DIMENSIONS	Height: 20.2" • Width: 10.9" • Depth: 8"
WEIGHT**	Device: 16.7 lbs • Battery: 1.35 lbs/each
OXYGEN FLOW RATE	Pulse Dose Settings 0.5-6 in 0.5 increments (8mL-96mL) 1-3 LPM in Continuous Mode in 0.5 LPM increments
BATTERY DURATION***	5.75 hours @ 2 pulse at 20 bpm • 2.5 hours @ 2 LPM continuous
SOUND LEVEL @ 2 PULSE	40 dBA
BATTERY RECHARGE TIME****	1.5 hours - 1 battery • 2.5 hours - 2 batteries
OXYGEN CONCENTRATION	95.6% - 87% at Sea Level at all settings
OXYGEN OUTPUT PRESSURE	5 psig +/- 1 psig (34 +/- 7 kPa)
AC POWER	Universal Input: 100-240 VAC • 50-60 Hz • 2.0 Amps Output: 24 Volts • 6.3 Amps
DC POWER	12 - 15 VDC • 10 Amps
OPERATING TEMPERATURE	50°F (10°C) to 104°F (40°C)
ALTITUDE OPERATING RANGE	0 - 13,123 Feet (0 - 4,000 Meters)
TRANSPORTATION/STORAGE TEMPERATURE	-4°F (-20°C) to 140°F (60°C)
TRANSPORTATION/STORAGE HUMIDITY	0 - 95% Non-Condensing
WARRANTY	3 Year Limited • 1 Year on Accessories/Battery



*O2 CONCEPTS RESERVES THE RIGHT TO CHANGE SPECIFICATIONS OF THIS PRODUCT AT ANYTIME WITHOUT NOTICE

**PRODUCTION WEIGHT DOES NOT INCLUDE CART KIT

***BATTERY DURATION SHOWN IS FOR EXTENDED BATTERY SYSTEM

****BATTERY RECHARGE TIME WILL INCREASE IF UNIT IS RUNNING WHILE CHARGING

ORDERING INFORMATION*

PART DESCRIPTION	PART NUMBER
SYSTEMS	
OXLIFE INDEPENDENCE POC STANDARD SYSTEM The Oxlife Independence POC Standard System includes oxygen concentrator, 1 battery, pull handle and wheels, AC and DC power supplies, accessory bag and user manual.	800-0001-1
OXLIFE INDEPENDENCE POC EXTENDED BATTERY SYSTEM The Oxlife Independence POC Extended Battery System includes oxygen concentrator, 2 batteries, pull handle and wheels, AC and DC power supplies, accessory bag and user manual.	800-0001-2
ACCESSORIES	
AC POWER CORD/POWER SUPPLY	800-1000
DC POWER CORD	800-1001
BATTERY - LITHIUM ION RECHARGEABLE	800-1002
ACCESSORY CARRY BAG	800-1003
BATTERY CHARGER - DESK TOP SINGLE BAY	800-1004
BATTERY CHARGER - MULTIPLE (6 BAY)	800-1005
INLET AIR FILTER (PACK OF 6)	800-1006
WHEELS (2)	800-1007
USER MANUAL	800-1008
PROTECTIVE COVER	800-1009
CART KIT (INCLUDES WHEELS AND PULL HANDLE)	800-1010
PULL HANDLE KIT	800-1011

*US FEDERAL LAW RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN

AUTHORIZED DEALER



Oxygen

Portable

USA

Your one stop shop for Portable Oxygen Equipment, Service and Accessories

Live Customer Support **1.800.339.1558**
9am to 6pm PST



PROUDLY MADE IN THE USA