

Portable Oxygen Concentrator Comparison Chart

	Airsep Focus	Respironics SimplyGO	Precision Medical Easy Pulse	Airsep Freestyle 5	Respironics EverGo	Lifechoice Activox	O2 Concepts Oxlife Independence	Devilbiss iGO
								
Rx Required	YES	YES	YES	YES	YES	YES	YES	YES
Flow Settings								
Continuous	Pulse Only	1-2 LPM	Pulse Only	Pulse Only	Pulse Only	Pulse Only	1-2-3 LPM	1-2-3 LPM
Pulse	2 LPM	1-6 LPM	1-5 LPM	1-5 LPM	1-6 LPM	1-3 LPM	1-6 LPM	1-6 LPM
Battery Duration								
Continuous								
1 LPM	NA	2.3 hours	NA	NA	NA	NA	3.5 hours	4 hours
2 LPM	NA	0.9 hours	NA	NA	NA	NA	2.5 hours	2.4 hours
3 LPM	NA	NA	NA	NA	NA	NA	1.5 hours	1.6 hours
Pulse								
1 LPM	NA	3.7 hours	4.5 hours	3 hours	6 hours	12.5 hours - Pro	5.4 hours	5.4 hours
2 LPM	3 hours	3.5 hours	3.2 hours	2 hours	4 hours	7.25 hours - Pro	4.7 hours	4.7 hours
3 LPM	NA	3.1 hours	2.5 hours	1.5 hours	3 hours	6.25 hours - Pro	4 hours	4 hours
4 LPM	NA	2.7 hours	1.8 hours	1.25 hours	3 hours	NA	3 hours	3.5 hours
5 LPM	NA	2.3 hours	1.5 hours	1 hour	2 hours	NA	2.25 hours	3.2 hours
6 LPM	NA	1.9 hours	NA	NA	2 hours	NA	1.75 hours	3.0 hours
O2 Concentration	90% to 95%	86% to 97%	87% to 95%	90% to 95%	87% to 93%	90% to 93%	87% to 95%	88 to 93%
Maximum Delivery	380 ml/min	2000 ml/min	780 ml/min	1000 ml/min	1050 ml/min	3000 ml/min	3000 ml /min	960 ml/min
Dimensions	6.4" x 4.8" x 2.5" H x W x D	11.5" x 10" x 6" H x W x D	10.1" x 6.5" x 4.5" H x W x D	10.7" x 6.6" x 4.4" H x W x D	8.5" x 6" x 12" H x W x D	7.87" x 9.05" x 4.38" H x W x D	20.2" x 10.9" x 8" H x W x D	15" x 11" x 8" H x W x D
Weight	2 lbs	10 lbs	6.8 lbs	6.2 lbs	10 lbs	3.9 lbs - Sport 4.3 lbs - Pro	16.7 lbs	20 lbs
Maximum Operational Altitude	up to 10,000 ft	up to 10,000 ft	up to 9,000 ft	up to 12,000 ft	up to 8,000 ft	up to 10,000 ft	up to 13,000 ft	up to 13,000 ft
Noise	41 dba	43 dba	44 dba	41 dba	43 dba	46 dba	40 dba	40 dba
Approximate Battery Recharge Time	4 hours	3.5 hours	4 hours	3 hours	3 hours	3 hours - Sport 4.5 hours - Pro	3.1 hours	4.5 hours
Car charger compatible	all settings	all settings	all settings	all settings	all settings	all settings	all settings	up to 2 LPM
Battery will charge in car	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No
Factory Warranty (Excluding Battery)	3 years	2 years	3 years	3 years	3 years	3 years	3 years	3 years limited
	Buy Now	Buy Now	Buy Now	Buy Now	Buy Now	Buy Now	Buy Now	Buy Now

All portable oxygen concentrators are FAA approved.
42 dba is the normal sound level of a human talking at a distance of 7.3 meters.