










# Portable Oxygen Concentrator Comparison Chart

ALL UNITS ARE FAA APPROVED



Your one stop shop for Portable Oxygen Equipment, Service and Accessories  
1.800.339.1558

	Airsep Focus	Precision Medical Easy Pulse	Airsep Freestyle 5	O2 Concepts Oxlife Independence	Sequal Eclipse 3	Devilbiss iGO	Respironics Evergo	Invacare XPO2
								
<b>Rx Required</b>	YES	YES	YES	YES	YES	YES	YES	YES
<b>Flow Settings</b>								
<b>Continuous</b>								
<b>Pulse</b>	Pulse Only 2 LPM	Pulse Only 1-5 LPM	Pulse Only 1-5 LPM	1-2-3 LPM 1-6 LPM	1-2-3 LPM 1-6 LPM	1-2-3 LPM 1-6 LPM	Pulse Only 1-6 LPM	Pulse Only 1-5 LPM
<b>Battery Duration</b>								
<b>Continuous</b>								
1 LPM	NA	NA	NA	3.5 hours	3.7 hours	4 hours	NA	NA
2 LPM	NA	NA	NA	2.5 hours	2 hours	2.4 hours	NA	NA
3 LPM	NA	NA	NA	1.5 hours	1.3 hours	1.6 hours	NA	NA
<b>Pulse</b>								
1 LPM	NA	4.5 hours	3 hours	5.4 hours	5.4 hours	5.4 hours	4 hours	7 hours
2 LPM	3 hours	3.2 hours	2 hours	4.7 hours	4.5 hours	4.7 hours	4 hours	5 hours
3 LPM	NA	2.5 hours	1.5 hours	4 hours	4 hours	4 hours	3 hours	4 hours
<b>O2 Concentration</b>	90% to 95%	87% to 95%	90% to 95%	87% to 95%	90% to 93%	88% to 94%	86% to 92%	87% to 95%
<b>Maximum Delivery</b>	380 ml/min	780 ml/min	1000 ml/min	3000 ml /min	3000 ml/min	960 ml/min	1050 ml/min	950 ml/min
<b>Dimensions</b>	6.4" x 4.8" x 2.5" H x W x D	10.1" x 6.5" x 4.5" H x W x D	10.7" x 6.6" x 4.4" H x W x D	20.2" x 10.9" x 8" H x W x D	19.3" x 12.3" x 7.1" H x W x D	15" x 11" x 8" H x W x D	8.5" x 6" x 12" H x W x D	10" x 7" x 4" H x W x D
<b>Weight</b>	2 lbs	6.8 lbs	6.2 lbs	16.7 lbs	20 lbs	20 lbs	10 lbs	6 lbs
<b>Maximum Operational Altitude</b>	up to 12,000 ft	up to 9,000 ft	up to 12,000 ft	up to 13,000 ft	up to 13,000 ft	up to 13,000 ft	up to 8,000 ft	up to 10,000 ft
<b>Car charger compatible</b>	all settings	all settings	all settings	all settings	up to 2 LPM	up to 2 LPM	all settings	all settings
<b>Battery will charge in car</b>	Yes	Yes	Yes	Yes	No	No	Yes	Yes
<b>Factory Warranty (Excluding Battery)</b>	3 years	3 years	3 years	3 years	3 years	3 years	3 years	3 years
	<b>Buy Now</b>	<b>Buy Now</b>	<b>Buy Now</b>	<b>Buy Now</b>	<b>Buy Now</b>	<b>Buy Now</b>	<b>Buy Now</b>	<b>Buy Now</b>